



i Roll Junior

Blood Tube Rotator



- AC synchronous motor for noiseless operation.
- Can accommodate upto 10 Blood collection tubes
- Circular motion gives gentle but effective mixing
- Fix speed of 30 RPM with an accuracy of ± 2 RPM
- Continuous / infinite mode of operation
- Robust construction with strong metal from design
- Maximum operational weight capacity of 1kg
- Compact saves bench space design
- Protection class of IP 21
- Easy to remove and install attachments for replacement and cleaning

Specification	iRoll Junior
Capacity (supplied with any one chosen attachments)	10 x 6 ml
	9 x 15 ml
	4 x 50 ml
Motor Type	AC Synchronous
Variable Speed	Fix Speed of 30 RPM
Run Time	Continuous / Infinite
Speed Accuracy	± 2 RPM
Removable Plate	Yes
Noise Level	Noiseless operation
Max Loading (Capacity)	1 kg
Dimension (W x D x H)	190 x 135 x 133 mm
Protection class	IP 21



* OEM request can be entertained. Specifications can be changed without notice for quality improvement.